



DENTAL

# Instrument Processing



# Midmark M9<sup>®</sup> + M11<sup>®</sup> Steam Sterilizers

Midmark sterilizers are designed to be safe, dependable and easy to use, featuring pre-programmed sterilization cycles, an intuitive LCD display with prompts to ensure quick, confident operation, and real-time safety alerts.

## SAFETY FEATURES

### Simple Control:


Once a pre-programmed cycle is selected, the unit is designed to sterilize automatically without any operator assistance.

### Smart Technology:

The sterilizer will automatically stop and alert the operator to take proper action if the door is not closed completely or the water level is low.

Model/Order No.	M9	M11
Overall length with plug	20.38" (51.8 cm)	23.8" (60.5 cm)
Overall width	15.3" (38.9 cm)	17.8" (45.2 cm)
Overall height with datalogger	17.24" (43.8 cm)	19.24" (48.9 cm)
Overall height with printer	16.1" (40.9 cm)	18.1" (50.0 cm)
Minimum countertop area	17.88" (45.4 cm) D x 15.3" (38.9 cm) W	21" (53.3 cm) D x 17.8" (45.2 cm) W
Chamber	9" diameter (22.9 cm) 15" depth (38.1 cm) 3.5 gal usable volume (13.4 liters)	11" diameter (28 cm) 18" depth (45.7 cm) 6.5 gal usable volume (24.6 liters)
Trays and pouch racks	Two large trays – each 7.3" W x 12" L x 0.8" D (18.6 cm x 30.5 cm x 2.2 cm) Two small trays – each 5.6" W x 12" L x 0.8" D (14.3 cm x 30.5 cm x 2.2 cm) One 5-slot pouch rack	Two large trays – each 9" W x 15" L x 1.1" D (22.9 cm x 38 cm x 2.9 cm) Two small trays – each 6.6" W x 15" L x 1.1" D (16.8 cm x 38 cm x 2.9 cm) One 6-slot pouch rack
Unit weight	73 lb (33.1 kg)	99 lb (44.9 kg)
Shipping weight	81 lb (36.7 kg)	131 lb (59.4 kg)
Water reservoir capacity	1.1 gal (4.1 liters)	1.4 gal (5.3 liters)
Electrical requirements	See midmark.com/IP for details	See midmark.com/IP for details

Steam sterilizers include a 1-year limited warranty.



Midmark M9 and M11 Steam Sterilizers include a 1-year limited warranty. Add a Sterilizer Data Logger at any time.

**01** Midmark sterilizers require no additional air or water filtration accessories or expensive daily Bowie-Dick tests. They cost only pennies per day to maintain compared to the dollars per day required by other units.

**02** Place your pouches in pouch racks to allow heat and steam to effectively circulate throughout the sterilizer chamber. Always use pouch racks on the bottom tray.

**03** All Midmark sterilizers are designed, manufactured and supported in the US.